

## Vaccum blackhead remover User Manual

"The vacuum blackhead remover" is a multifunctional beauty device designed to extract blackheads from the face, by a principle of operation "unique suction system" exfoliating, firming and gently massaging the skin of the face. It results a better blood and lymphatic circulation, to make it more compact, thus finding a healthier and brighter skin.

Safe and non-toxic, made of eco-friendly ABS material, the blackhead machine is reliable with **CE Fcc RoHS** certification, guaranteeing your safety when using it.

### **PRECAUTION FOR USE:**

#### **IMPORTANT:**

A touch test is required before any use of the multifunctional beauty device. Test it on your entire forearm, above and below, for at least 15 minutes.

If you do not notice any reaction on your skin, you can use it.

For dry or damaged skin, it is better to use it the bare minimum, once a week. Stimulates the blood and the lymphatic system for a healthy and glowing skin. Used for dull and tired skin.

Do not stay more than 3 seconds on the surface to be treated with the tips (tips), small hickeys or hematomas should occur. In the event that the use of the machine causes small bruises or small purplish marks, let your skin rest, and start again the next day. For fine and sensitive skin, it is preferable to use tip No. 2: small circular tip, once a week.

After each session, **BE CAREFUL**, do not apply exfoliating or astringent products or those used to remove dead skin during the next 3 days following your care. The blackhead machine (multifunctional beauty device) is in no way a miracle product, it is a multifunctional beauty device that complements your weekly care, do not overuse it.

Using it once or twice a week is enough.

Do not leave within reach of young children.

The use of the blackhead machine (multifunctional beauty device) is not appropriate for children under 13 years old.

This device is for strictly personal use only, do not lend it. What is appropriate for you may not be convenient for others, even a member of your family.

In case of doubt, it is asked to take advice from your usual specialist.

## **USE:**

- 1- Charge your device, for first use charge it with the USB cable for about 6 hours and wait until the LED is green.
- 2- Take the multifunctional beauty device and fit the silicone ring into the notch provided for this purpose at the top of your machine.
- 3- To properly dilate the pores of the skin and facilitate the extraction of blackheads, apply a warm towel or a warm glove on your face for 2 minutes, even more, depending on your skin type.
- 4- Wipe your face (no water in the device).
- 5- Take the device and put the first tip big size circular.
- 6- Turn it on and slide the tip up and down on 1 to 2 cms each time, but do not stay on the same surface for more than 3 seconds.
- 7- Make all areas of your face (except the eyelids and the underside of the eyes). Do not press too hard, redness may appear on the skin.
- 8- Continue with the other tips according to your expectations (see explanation of the tips).

**Once you have finished using your multifunctional beauty device, you can apply an extra-gentle mask or cream to complete your care.**

Use the device once or twice a week, (except for very dry, dry, damaged, fine and sensitive skin). After a while, you will notice an increase of the skin elasticity, a better blood circulation in the face and that your skin is regaining radiance. For more results, make a facial sauna to dilate the pores of the skin before your session with the blackhead machine. Do not fully empty the device to save battery and keep your machine for longer.

## **CLEANING THE DEVICE AND THE TIPS:**

Due to the extraction of the blackheads with the tips, bacteria can agglutinate therein. It is strongly advised to sterilize the tips in boiling water or to buy disinfecting wipes for medical or aesthetic equipment.

## **REPLACEMENT FILTERS**

Remove the suction head then the damaged filter. Do not immerse the device.

## INFORMATION:

The ACDZ company headquartered **183 Chemin des Négadoux 83140 Six - Fours-les – Plages** is not and will not be held responsible for any other use or non-compliance with the precaution of use mentioned in this booklet and asks users to respect them.

The “Blackhead machine” (multifunctional beauty device) is 1 year warranty. In case of dispute you should refer to the Terms of Sale of NAYCIA.FR You can contact the after-sales service at: **contact@naycia.fr**

## CHARACTERISTICS:

- Made of environmentally friendly BS materials
- Four different non-abrasive heads / tips (made of a microcrystalline diamond tip)
- Lithium polymer battery battery type 450Ah – 3.7V
- Weight: 127.4 gr - Length: 18.8 X4.5 X4 cm

## CONTAINS

- 1 Vacuum cleaner**
- 1 USB cable**
- 4 Non-abrasive heads/probes**
- 4 replacement filters**
- (Avoid impurities to enter the device)**
- 2 silicon rings (For vacuum heads stabilization)**

### 1- Large circular size tip :

Designed to effectively remove blackheads by strong suction. Suitable for all areas of the face except the eyelids and the underside of the eyes.

### 2 - Circular small size tip:

It is best to use this one for fine and sensitive skins. Suitable for all areas of the face except the eyelids and the underside of the eyes. Softer and more delicate, it is designed to unclog blackheads when the large tip is unable to do it (pores too tight).

### 3- Elliptical tip:

Designed to make the skin more elastic and compact.

Used to remove small wrinkles on the face.

Sessions with this elliptical tip restore the elasticity of the skin. Suitable for all areas of the face except the eyelids and the underside of the eyes.

#### **4- Micro-Crystal tip:**

Used to exfoliate keratin and remove dead skin from the face. Stimulates the blood and the lymphatic system, helps to regain a healthy and luminous skin. Used for dull and tired skins. Suitable for all areas of the face except the eyelids and the underside of the eyes.

#### **SELECTIVE SORTING OF THE BATTERY**

- The built-in rechargeable battery contains substances that can harm the environment.
- Before disposing off this unit, the rechargeable battery must be removed.
- The device must be disconnected from the power supply network when removing the battery.
- The battery must be disposed of safely.
- Do not dispose off the battery or appliance with household wastes. To preserve the environment, dispose off the device and the battery in accordance with the regulations in force.
- Deposit them in the collection bins provided for this purpose or bring them back to your store. Never disassemble, crush or pierce the battery.
- Do not throw it into the fire.
- Do not expose it to high temperatures or direct sunlight.
- Do not incinerate it.